

Phospho-JUND (S255) Antibody



PACO05955

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IP

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IP:2-5ug

Protein Background:

Gene ID:

JUND

Uniprot

P17535

Synonyms:

JUND; Transcription factor jun-D

Immunogen:

Synthesized peptide derived from human Jun D around the phosphorylation site of S255.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A