

NDUFS7 Antibody



PACO05957

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:**Gene ID:**

NDUFS7

Uniprot

O75251

Synonyms:

NDUFS7; NADH dehydrogenase [ubiquinone] iron-sulfur protein 7; mitochondrial; Complex I-20kD; CI-20kD; NADH-ubiquinone oxidoreductase 20 kDa subunit; PSST subunit

Immunogen:

Synthesized peptide derived from the C-terminal region of human NDUFS7.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A