

FAM111B Antibody



PACO05958

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:

Gene ID:

FAM111B

Uniprot

Q6SJ93

Synonyms:

FAM111B; CANP; Protein FAM111B; Cancer-associated nucleoprotein

Immunogen:

Synthesized peptide derived from the Internal region of human FAM111B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A