

KCND2 Antibody



PACO05966

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300

Protein Background:**Gene ID:**

KCND2

Uniprot

Q9NZV8

Synonyms:

KCND2; KIAA1044; Potassium voltage-gated channel subfamily D member 2; Voltage-gated potassium channel subunit Kv4.2

Immunogen:

Synthesized peptide derived from human Kv4.2 around the non-phosphorylation site of S616.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A