

KCNF1 Antibody



PACO05967

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:**Gene ID:**

KCNF1

Uniprot

Q9H3M0

Synonyms:

KCNF1; Potassium voltage-gated channel subfamily F member 1; Voltage-gated potassium channel subunit Kv5.1; KH1

Immunogen:

Synthesized peptide derived from the Internal region of human KCNF1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A