

# MAST4 Antibody



PACO05968

---

## Product Information

<b>Size:</b> 50ug	<b>Protein Background:</b>
<b>Reactivity:</b> Human, Mouse	<b>Gene ID:</b> MAST4
<b>Source:</b> Rabbit	<b>Uniprot</b> O15021
<b>Isotype:</b> IgG	<b>Synonyms:</b> MAST4; KIAA0303; Microtubule-associated serine/threonine-protein kinase 4
<b>Applications:</b> ELISA, WB, IHC, IF	<b>Immunogen:</b> Synthesized peptide derived from the N-terminal region of human MAST4.
<b>Recommended dilutions:</b> ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000	<b>Storage:</b> Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A