

Phospho-KCNJ1 (S44) Antibody



PACO05969

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

KCNJ1

Uniprot

P48048

Synonyms:

KCNJ1; ROMK1; ATP-sensitive inward rectifier potassium channel 1; ATP-regulated potassium channel ROM-K; Inward rectifier K(+) channel Kir1.1; Potassium channel; inwardly rectifying subfamily J member 1

Immunogen:

Synthesized peptide derived from human ROM-K around the phosphorylation site of S44.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A