

OR5B3 Antibody



PACO05979

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:10000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR5B3

Uniprot

Q8NH48

Synonyms:

OR5B3; OR5B13; Olfactory receptor 5B3; Olfactory receptor 5B13; Olfactory receptor OR11-239

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 5B3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A