

Phospho-HDAC3 (S424) Antibody



PACO05998

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

HDAC3

Uniprot

O15379

Synonyms:

HDAC3; Histone deacetylase 3; HD3; RPD3-2; SMAP45

Immunogen:

Synthesized peptide derived from human HDAC3 around the phosphorylation site of S424.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A