

OR51B5 Antibody



PACO06006

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR51B5

Uniprot

Q9H339

Synonyms:

OR51B5; Olfactory receptor 51B5; Odorant receptor HOR5'beta5; Olfactory receptor OR11-37

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 51B5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A