

Acetyl-FOXA1 (K264/253/211) Antibody



PACO06058

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

FOXA1

Uniprot

P55317

Synonyms:

FOXA1; HNF3A; TCF3A; Hepatocyte nuclear factor 3-alpha; HNF-3-alpha; HNF-3A; Forkhead box protein A1; Transcription factor 3A; TCF-3A; FOXA2; HNF3B; TCF3B; Hepatocyte nuclear factor 3-beta; HNF-3-beta; HNF-3B;


Immunogen:

Synthesized peptide derived from the human HNF-3 α / β / γ ; around the acetylation site of K264/253/211.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Acetyl-HNF-3alpha/beta/gamma(K264/253/211) Polyclonal Antibody.



Western Blot analysis of JK cells using Acetyl-HNF-3alpha/β/gamma(K264/253/211) Polyclonal Antibody.