

Acetyl-Histone H2B (K126) Antibody



PACO06059

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

HIST1H2BA

Uniprot

Q96A08

Synonyms:

HIST1H2BA; TSH2B; Histone H2B type 1-A; Histone H2B, testis; Testis-specific histone H2B; HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B. f; H2B/f; HIST1H2BC; H2BFL; HIST1H2BE; H2BFH; HIST1H2BF;


Immunogen:

Synthesized peptide derived from the human Histone H2B around the acetylation site of K126.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Acetyl-Histone H2B (K126) Polyclonal Antibody.



Western Blot analysis of A549 cells using Acetyl-Histone H2B (K126) Polyclonal Antibody.