

Acetyl-HSP90AA1 (K292/284) Antibody



PACO06061

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

HSP90AA1

Uniprot

P07900

Synonyms:

HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38; HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84


Immunogen:

Synthesized peptide derived from the human HSP 90 around the acetylation site of K292/284.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using Acetyl-HSP 90 (K292/284) Polyclonal Antibody.



Western Blot analysis of K562 cells using Acetyl-HSP 90 (K292/284) Polyclonal Antibody.