

Acetyl-TP73 (K327) Antibody



PACO06068

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

TP73

Uniprot

O15350

Synonyms:

TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein


Immunogen:

Synthesized peptide derived from the human p73 around the acetylation site of K327.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Acetyl-p73 (K327) Polyclonal Antibody.



Western Blot analysis of JK cells using Acetyl-p73 (K327) Polyclonal Antibody.