

# Acetyl-UBA52/RPS27A/UBB/UBC (K29) Antibody



PACO06079

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000

**Protein Background:****Gene ID:**

UBA52/RPS27A/UBB/UBC

**Uniprot**

P62987/P62979/P0CG47/P0CG48

**Synonyms:**

UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A; UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2; Ubiquitin-60S ribosomal protein L40; CEP52; Ubiquitin A-52 residue ribosomal protein fusion pr

**Immunogen:**


Synthesized peptide derived from the human Ub around the acetylation site of K29.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of HepG2 cells using Acetyl-Ub (K29) Polyclonal Antibody.



Western Blot analysis of HepG2 cells using Acetyl-Ub (K29) Polyclonal Antibody.