

Acetyl-UBA52/RPS27A/UBB/UBC (K33) Antibody



PACO06080

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

UBA52/RPS27A/UBB/UBC

Uniprot

P62987/P62979/P0CG47/P0CG48

Synonyms:

UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A; UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2; Ubiquitin-60S ribosomal protein L40; CEP52; Ubiquitin A-52 residue ribosomal protein fusion pr


Immunogen:

Synthesized peptide derived from the human Ub around the acetylation site of K33.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Acetyl-Ub (K33) Polyclonal Antibody.



Western Blot analysis of HeLa cells using Acetyl-Ub (K33) Polyclonal Antibody.