
Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Rabbit, Chicken,
Monkey, Sheep, X

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

WB:1:3000, IHC:1:200

Protein Background:**Gene ID:**

ACTB

Uniprot

P60709

Synonyms:

ACTB; Actin, cytoplasmic 1; beta-actin


Immunogen:

Recombinant Protein.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of HepG2 (1), Rat liver (2), Mouse kidney (3), Rabbit testis (4), Sheep lung (5), 293T (6), diluted at 1:5000.



IHC staining of Human ovary tissue, diluted at 1:200.