

Phospho-JUN/JUND (S73/100) Antibody



PACO06094

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

JUN/JUND

Uniprot

P05412/P17535

Synonyms:

JUN; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun
avian sarcoma virus 17 oncogene homolog; p39; JUND; Transcription factor jun-D


Immunogen:

Synthesized peptide derived from human AP-1 around the phosphorylation site of S73.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of VEC A549 cells using Phospho-AP-1/Jun D (S73/100) Polyclonal Antibody.



Western Blot analysis of HuvEc cells using Phospho-AP-1/Jun D (S73/100) Polyclonal Antibody.