

Phospho-Histone H2A. X (S139) Antibody



PACO06127

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

H2AFX

Uniprot

P16104

Synonyms:

H2AFX; H2AX; Histone H2A. x; H2a/x


Immunogen:

Synthesized peptide derived from human Histone H2A. X around the phosphorylation site of S139.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Phospho-Histone H2A. X (S139) Polyclonal Antibody.



Western Blot analysis of 4T1 VEC 293 cells using Phospho-Histone H2A. X (S139) Polyclonal Antibody.