

Phospho-MEF2A (T319) Antibody



PACO06139

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

MEF2A

Uniprot

Q02078

Synonyms:

MEF2A; MEF2; Myocyte-specific enhancer factor 2A; Serum response factor-like protein 1


Immunogen:

Synthesized peptide derived from human MEF-2 around the phosphorylation site of T319.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of VEC 293 cells using Phospho-MEF-2 (T319) Polyclonal Antibody.



Western Blot analysis of VEC cells using Phospho-MEF-2 (T319) Polyclonal Antibody.