

# Phospho-TP73 (Y99) Antibody



PACO06148

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

TP73

**Uniprot**

O15350

**Synonyms:**

TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein

**Immunogen:**


Synthesized peptide derived from human p73 around the phosphorylation site of Y99.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of NIH-3T3 cells using Phospho-p73 (Y99) Polyclonal Antibody.



Western Blot analysis of HeLa cells using Phospho-p73 (Y99) Polyclonal Antibody.