

# Phospho-HDAC2 (S394) Antibody



PACO06156

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat, Monkey

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

HDAC2

**Uniprot**

Q92769

**Synonyms:**

HDAC2; Histone deacetylase 2; HD2

**Immunogen:**


Synthesized peptide derived from human HDAC2 around the phosphorylation site of S394.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of NIH-3T3 cells using Phospho-HDAC2 (S394) Polyclonal Antibody.



Western Blot analysis of 3T3 cells using Phospho-HDAC2 (S394) Polyclonal Antibody.