

# ABL1/ABL2 Antibody



PACO06221

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:****Gene ID:**

ABL1/ABL2

**Uniprot**

P00519/P42684

**Synonyms:**

ABL1; ABL; JTK7; Tyrosine-protein kinase ABL1; Abelson murine leukemia viral oncogene homolog 1; Abelson tyrosine-protein kinase 1; Proto-oncogene c-Abl; p150; ABL2; ABLL; ARG; Abelson tyrosine-protein kinase 2; Abelson murine leukemia vira

**Immunogen:**


Synthesized peptide derived from human Abl1/2 around the non-phosphorylation site of Y393/439.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of K562 cells using Abl1/2 Polyclonal Antibody.



Western Blot analysis of K562 cells using Abl1/2 Polyclonal Antibody.