

ACAT1 Antibody



PACO06222

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

ACAT1

Uniprot

P24752

Synonyms:

ACAT1; ACAT; MAT; Acetyl-CoA acetyltransferase; mitochondrial; Acetoacetyl-CoA thiolase; T2


Immunogen:

Synthesized peptide derived from the Internal region of human ACAT-1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc MCF7 HELA cells using ACAT-1 Polyclonal Antibody.



Western Blot analysis of A549 cells using ACAT-1 Polyclonal Antibody.