

ACTN1/ACTN2/ACTN3/ACTN4 Antibody



PACO06230

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

ACTN1/ACTN2/ACTN3/ACTN4

Uniprot

P12814/P35609/Q08043/O43707

Synonyms:

ACTN1; alpha-actinin-1; alpha-actinin cytoskeletal isoform; F-actin cross-linking protein; Non-muscle alpha-actinin-1; ACTN2; alpha-actinin-2; alpha-actinin skeletal muscle isoform 2; F-actin cross-linking protein; ACTN3; alpha-actinin-3; A


Immunogen:

Synthesized peptide derived from the N-terminal region of human Actinin- α 1/2/3/4.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat 293 HELA A549 cells using Actinin- α 1/2/3/4 Polyclonal Antibody.



Western Blot analysis of 293 cells using Actinin- α 1/2/3/4 Polyclonal Antibody.