

AIFM1 Antibody



PACO06239

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

AIFM1

Uniprot

O95831

Synonyms:

AIFM1; AIF; PDCD8; Apoptosis-inducing factor 1; mitochondrial; Programmed cell death protein 8


Immunogen:

Synthesized peptide derived from the N-terminal region of human AIF-M1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc HELA cells using AIF-M1 Polyclonal Antibody.



Western Blot analysis of Jurkat cells using AIF-M1 Polyclonal Antibody.