

# JUN/JUND Antibody



PACO06253

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

JUN/JUND

**Uniprot**

P05412/P17535

**Synonyms:**

JUN; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun  
avian sarcoma virus 17 oncogene homolog; p39; JUND; Transcription factor jun-D

**Immunogen:**


Synthesized peptide derived from human AP-1 around the non-phosphorylation site of  
S73.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of NIH-3T3 cells using AP-1/Jun D Polyclonal Antibody.



Western Blot analysis of HuvEc cells using AP-1/Jun D Polyclonal Antibody.