

APAF1 Antibody



PACO06254

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

APAF1

Uniprot

O14727-6

Synonyms:

APAF1-interacting protein; APIP; APIP2; CGI-29 protein; MMRP19; dJ179L10.2; likely ortholog of mouse monocyte macrophage 19


Immunogen:

Synthesized peptide derived from the C-terminal region of human Apaf-1-ALT.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Apaf-1-ALT Polyclonal Antibody.



Western Blot analysis of COLO205 cells using Apaf-1-ALT Polyclonal Antibody.