

APC Antibody



PACO06255

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

APC

Uniprot

P25054

Synonyms:

APC; DP2.5; Adenomatous polyposis coli protein; Protein APC; Deleted in polyposis 2.5


Immunogen:

Synthesized peptide derived from the C-terminal region of human APC.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-MUSCLE cells using APC Polyclonal Antibody.



Western Blot analysis of COLO205 cells using APC Polyclonal Antibody.