

ACTR3 Antibody



PACO06264

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

ACTR3

Uniprot

P61158

Synonyms:

ACTR3; ARP3; Actin-related protein 3; Actin-like protein 3


Immunogen:

Synthesized peptide derived from the C-terminal region of human Arp3.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using Arp3 Polyclonal Antibody.



Western Blot analysis of K562 cells using Arp3 Polyclonal Antibody.