
Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

ACLY

Uniprot

P53396

Synonyms:

ACLY; ATP-citrate synthase; ATP-citrate; pro-S-)-lyase; ACL; Citrate cleavage enzyme


Immunogen:

Synthesized peptide derived from human ATP-citrate synthase around the non-phosphorylation site of S455.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293T cells using ATP-citrate synthase Polyclonal Antibody.



Western Blot analysis of COS7 cells using ATP-citrate synthase Polyclonal Antibody.