

MCM2 Antibody



PACO06288

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

MCM2

Uniprot

P49736

Synonyms:

MCM2; BM28; CCNL1; CDCL1; KIAA0030; DNA replication licensing factor MCM2;
Minichromosome maintenance protein 2 homolog; Nuclear protein BM28


Immunogen:

Synthesized peptide derived from the N-terminal region of human BM28.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF7 PC-3 KB HELA cells using BM28 Polyclonal Antibody.



Western Blot analysis of 293 cells using BM28 Polyclonal Antibody.