

BMP7 Antibody



PACO06292

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

BMP7

Uniprot

P18075

Synonyms:

BMP7; OP1; Bone morphogenetic protein 7; BMP-7; Osteogenic protein 1; OP-1;
Eptoterminal alpha


Immunogen:

Synthesized peptide derived from the Internal region of human BMP-7.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-MUSCLE cells using BMP-7 Polyclonal Antibody.



Western Blot analysis of rat-muscle cells using BMP-7 Polyclonal Antibody.