

SLC4A11 Antibody



PACO06303

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

SLC4A11

Uniprot

Q8NBS3

Synonyms:

SLC4A11; BTR1; Sodium bicarbonate transporter-like protein 11; Bicarbonate transporter-related protein 1; Sodium borate cotransporter 1; NaBC1; Solute carrier family 4 member 11


Immunogen:

Synthesized peptide derived from the Internal region of human BTR1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using BTR1 Polyclonal Antibody.



Western Blot analysis of K562 cells using BTR1 Polyclonal Antibody.