

BUB3 Antibody



PACO06304

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

BUB3

Uniprot

O43684

Synonyms:

BUB3; Mitotic checkpoint protein BUB3


Immunogen:

Synthesized peptide derived from the C-terminal region of human Bub3.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Bub3 Polyclonal Antibody.



Western Blot analysis of 293 cells using Bub3 Polyclonal Antibody.