

CDC37 Antibody



PACO06337

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

CDC37

Uniprot

Q16543

Synonyms:

CDC37; CDC37A; Hsp90 co-chaperone Cdc37; Hsp90 chaperone protein kinase-targeting subunit; p50Cdc37


Immunogen:

Synthesized peptide derived from human Cdc37 around the non-phosphorylation site of S13.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using Cdc37 Polyclonal Antibody.



Western Blot analysis of MCF7 cells using Cdc37 Polyclonal Antibody.