

COL5A3 Antibody



PACO06368

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

COL5A3

Uniprot

P25940

Synonyms:

COL5A3; Collagen alpha-3(V) chain


Immunogen:

Synthesized peptide derived from the N-terminal region of human COL5A3.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using COL5A3 Polyclonal Antibody.



Western Blot analysis of K562 cells using COL5A3 Polyclonal Antibody.