

PACO06374

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

CKM

Uniprot

P06732

Synonyms:

CKM; CKMM; Creatine kinase M-type; Creatine kinase M chain; M-CK


Immunogen:

Synthesized peptide derived from the N-terminal region of human Creatine Kinase M.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Creatine Kinase M Polyclonal Antibody.



Western Blot analysis of Jurkat cells using Creatine Kinase M Polyclonal Antibody.