

CREB3L2 Antibody



PACO06376

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

CREB3L2

Uniprot

Q70SY1

Synonyms:

CREB3L2; BBF2H7; Cyclic AMP-responsive element-binding protein 3-like protein 2; cAMP-responsive element-binding protein 3-like protein 2; BBF2 human homolog on chromosome 7


Immunogen:

Synthesized peptide derived from the Internal region of human CREB3L2.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using CREB3L2 Polyclonal Antibody.



Western Blot analysis of HeLa cells using CREB3L2 Polyclonal Antibody.