

CYP19A1 Antibody



PACO06396

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

CYP19A1

Uniprot

P11511

Synonyms:

CYP19A1; ARO1; CYAR; CYP19; Cytochrome P450 19A1; Aromatase; CYPXIX;
Cytochrome P-450AROM; Estrogen synthase


Immunogen:

Synthesized peptide derived from the Internal region of human CYP19A1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using CYP19A1 Polyclonal Antibody.



Western Blot analysis of 293 cells using CYP19A1 Polyclonal Antibody.