

CYP2C19 Antibody



PACO06400

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

CYP2C19

Uniprot

P33261

Synonyms:

CYP2C19; Cytochrome P450 2C19; (R)-limonene 6-monooxygenase; (S)-limonene 6-monooxygenase; (S)-limonene 7-monooxygenase; CYP11C17; CYP11C19; Cytochrome P450-11A; Cytochrome P450-254C; Mephenytoin 4-hydroxylase


Immunogen:

Synthesized peptide derived from the Internal region of human CYP2C19.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA 293T cells using CYP2C19 Polyclonal Antibody.



Western Blot analysis of MOUSE-BRAIN cells using CYP2C19 Polyclonal Antibody.