

CYC1 Antibody



PACO06408

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

CYC1

Uniprot

P08574

Synonyms:

CYC1; Cytochrome c1; heme protein, mitochondrial; Complex III subunit 4; Complex III subunit IV; Cytochrome b-c1 complex subunit 4; Ubiquinol-cytochrome-c reductase complex cytochrome c1 subunit; Cytochrome c-1


Immunogen:

Synthesized peptide derived from the Internal region of human Cytochrome c1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-MUCLE cells using Cytochrome c1 Polyclonal Antibody.



Western Blot analysis of rat-muscle cells using Cytochrome c1 Polyclonal Antibody.