

POLB Antibody



PACO06431

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

POLB

Uniprot

P06746

Synonyms:

POLB; DNA polymerase beta


Immunogen:

Synthesized peptide derived from the C-terminal region of human DNA pol β .


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using DNA POL eta Polyclonal Antibody.



Western Blot analysis of K562 cells using DNA pol beta Polyclonal Antibody.