

ECI2 Antibody



PACO06438

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

ECI2

Uniprot

O75521

Synonyms:

ECI2; DRS1; HCA88; PECI; Enoyl-CoA delta isomerase 2; mitochondrial; DRS-1; Delta(3), delta(2)-enoyl-CoA isomerase; D3, D2-enoyl-CoA isomerase; Diazepam-binding inhibitor-related protein 1; DBI-related protein 1; Dodecenoyl-CoA isomerase; Hep


Immunogen:

Synthesized peptide derived from the Internal region of human DRS-1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using DRS-1 Polyclonal Antibody.



Western Blot analysis of 293 cells using DRS-1 Polyclonal Antibody.