

SLC1A1 Antibody



PACO06443

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

SLC1A1

Uniprot

P43005

Synonyms:

SLC1A1; EAAC1; EAAT3; Excitatory amino acid, transporter 3; Excitatory amino-acid, carrier 1; Neuronal and epithelial glutamate transporter; Sodium-dependent glutamate/aspartate transporter 3; Solute carrier family 1 member 1


Immunogen:

Synthesized peptide derived from the Internal region of human EAAT3.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using EAAT3 Polyclonal Antibody.



Western Blot analysis of 293 cells using EAAT3 Polyclonal Antibody.