

EIF2B3 Antibody



PACO06454

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:**Gene ID:**

EIF2B3

Uniprot

Q9NR50

Synonyms:

EIF2B3; Translation initiation factor eIF-2B subunit gamma; eIF-2B GDP-GTP exchange factor subunit gamma


Immunogen:

Synthesized peptide derived from the Internal region of human eIF2B γ .


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using eIF2B gamma Polyclonal Antibody.



Western Blot analysis of K562 cells using eIF2B gamma Polyclonal Antibody.