

PTGER4 Antibody



PACO06462

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

PTGER4

Uniprot

P35408

Synonyms:

PTGER4; PTGER2; Prostaglandin E2 receptor EP4 subtype; PGE receptor EP4 subtype; PGE2 receptor EP4 subtype; Prostanoid EP4 receptor


Immunogen:

Synthesized peptide derived from the Internal region of human EP4.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using EP4 Polyclonal Antibody.



Western Blot analysis of COLO205 cells using EP4 Polyclonal Antibody.