

FZD3 Antibody



PACO06490

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

FZD3

Uniprot

Q9NPG1

Synonyms:

FZD3; Frizzled-3; Fz-3; hFz3


Immunogen:

Synthesized peptide derived from the Internal region of human Frizzled-3.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Frizzled-3 Polyclonal Antibody.



Western Blot analysis of RAW264.7 cells using Frizzled-3 Polyclonal Antibody.