

FSHR Antibody



PACO06493

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

FSHR

Uniprot

P23945

Synonyms:

FSHR; LGR1; Follicle-stimulating hormone receptor; FSH-R; Follitropin receptor


Immunogen:

Synthesized peptide derived from the Internal region of human FSHR.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using FSHR Polyclonal Antibody.



Western Blot analysis of KB cells using FSHR Polyclonal Antibody.