

GABRA6 Antibody



PACO06499

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

GABRA6

Uniprot

Q16445

Synonyms:

GABRA6; gamma-aminobutyric acid, receptor subunit alpha-6; GABA(A) receptor subunit alpha-6


Immunogen:

Synthesized peptide derived from the Internal region of human GABAA R α 6.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MOUSE-BRAIN cells using GABAA R alpha 6 Polyclonal Antibody.



Western Blot analysis of MOUSE-BRAIN cells using GABAA Ralpha 6 Polyclonal Antibody.